

Product Datasheet

Bovine Serum Albumin Antibody (D1) - Chimeric - Azide and BSA Free NBP3-20114-0.2mg

Unit Size: 0.2 mg

Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-20114

Updated 3/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-20114



NBP3-20114-0.2mg

Bovine Serum Albumin Antibody (D1) - Chimeric - Azide and BSA Free

Product Information

Unit Size	0.2 mg
Concentration	1 mg/ml
Storage	Store at 4C for up to 3 months. For longer storage, aliquot and store at -20C.
Clonality	Monoclonal
Clone	D1
Preservative	0.02% Proclin 300
Isotype	IgG Kappa
Purity	Protein A purified
Buffer	PBS

Product Description


Host	Rabbit
Gene ID	280717
Gene Symbol	ALB
Species	Bovine
Immunogen	D1 was prepared by isolating anti-p24 scFvs from a phage display library made from a mouse immunized with lysate from HIV-1 infected H9 cells.

Product Application Details

Applications	ELISA, Immunoprecipitation
Recommended Dilutions	ELISA, Immunoprecipitation
Application Notes	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.

Images

ELISA: Bovine Serum Albumin Antibody (D1) - Chimeric - Azide and BSA Free [NBP3-20114] - Binding curves of the rabbit IgG chimeric version of the anti-BSA antibody D1 (NBP3-20114; blue line) and isotype control (NBP2-52638, anti-fluorescein; orange line) to an ELISA plate coated with BSA at a concentration of 5 ug/ml. A 3-fold serial dilution from 90 ug/ml was tested in the experiment. For signal detection, a 1:4000 dilution of an HRP-conjugated anti-rabbit IgG antibody was used.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-20114-0.2mg

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP1-96904	Bovine Serum Albumin Native Protein
NBP1-96755	Bovine Serum Albumin Native Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-20114

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

