

# Product Datasheet

## **TDRD3 Antibody (RAB-C419) [DyLight 755] - Chimeric NBP3-20149IR**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-20149IR](http://www.novusbio.com/NBP3-20149IR)

Updated 11/12/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-20149IR](http://www.novusbio.com/reviews/destination/NBP3-20149IR)



**NBP3-20149IR**

TDRD3 Antibody (RAB-C419) [DyLight 755] - Chimeric

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | RAB-C419                                                                                    |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG Kappa                                                                                   |
| Conjugate           | DyLight 755                                                                                 |
| Purity              | Protein A purified                                                                          |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                                            |
| Gene ID             | 81550                                                                                                                             |
| Gene Symbol         | TDRD3                                                                                                                             |
| Species             | Human                                                                                                                             |
| Immunogen           | This antibody was obtained by recombinant antibody (rAb) phage display recognizing TDRD3 protein under non-denaturing conditions. |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                 |

| Product Application Details |                                                                                                                                                           |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Mass Spectrometry |
| Recommended Dilutions       | Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Mass Spectrometry, Chromatin Immunoprecipitation (ChIP) |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-20149IR**

---

|               |                                   |
|---------------|-----------------------------------|
| NBP2-68850PEP | TDRD3 Recombinant Protein Antigen |
| NBP2-06186    | TDRD3 Overexpression Lysate       |
| NB100-1936    | SMN Antibody (2B1)                |
| NB100-56374   | Carm1 Antibody                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-20149IR](http://www.novusbio.com/reviews/submit/NBP3-20149IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

