

# Product Datasheet

## Dystrophin Antibody (DMD/6270) [DyLight 405] NBP3-20751V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-20751V](http://www.novusbio.com/NBP3-20751V)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-20751V](http://www.novusbio.com/reviews/destination/NBP3-20751V)



**NBP3-20751V**

Dystrophin Antibody (DMD/6270) [DyLight 405]

| Product Information         |                                                                                                |
|-----------------------------|------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                         |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.    |
| Storage                     | Store at 4C in the dark.                                                                       |
| Clonality                   | Monoclonal                                                                                     |
| Clone                       | DMD/6270                                                                                       |
| Preservative                | 0.05% Sodium Azide                                                                             |
| Isotype                     | IgG2a Kappa                                                                                    |
| Conjugate                   | DyLight 405                                                                                    |
| Purity                      | Protein A or G purified                                                                        |
| Buffer                      | 50mM Sodium Borate                                                                             |
| Product Description         |                                                                                                |
| Host                        | Mouse                                                                                          |
| Gene ID                     | 1756                                                                                           |
| Gene Symbol                 | DMD                                                                                            |
| Species                     | Human, Canine                                                                                  |
| Reactivity Notes            | Does not react with mouse, xenopus or chicken                                                  |
| Immunogen                   | Recombinant fragment (around aa1700-2300) of human Dystrophin (exact sequence is proprietary). |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.              |
| Product Application Details |                                                                                                |
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin            |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin            |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                         |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-20751V**

---

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| NBP1-96981V        | Mouse IgG2a Kappa Isotype Control (M2AK) [DyLight 405] |
| H00001756-P01-10ug | Recombinant Human Dystrophin GST (N-Term) Protein      |
| 291-G1-200         | IGF-I/IGF-1 [Unconjugated]                             |
| NBP2-10564         | Dystrophin Overexpression Lysate                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-20751V](http://www.novusbio.com/reviews/submit/NBP3-20751V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

