

# Product Datasheet

## SNX2 Antibody NBP3-21431-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-21431](http://www.novusbio.com/NBP3-21431)

Updated 5/8/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-21431](http://www.novusbio.com/reviews/destination/NBP3-21431)




**NBP3-21431-100ul****SNX2 Antibody**

| <b>Product Information</b>         |                                                                                                                                                                                         |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 100 ul                                                                                                                                                                                  |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                             |
| <b>Storage</b>                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                                  |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                              |
| <b>Preservative</b>                | 0.03% Proclin 300                                                                                                                                                                       |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                     |
| <b>Purity</b>                      | Antigen and protein A Affinity-purified                                                                                                                                                 |
| <b>Buffer</b>                      | PBS, pH 7.0                                                                                                                                                                             |
| <b>Product Description</b>         |                                                                                                                                                                                         |
| <b>Host</b>                        | Rabbit                                                                                                                                                                                  |
| <b>Gene ID</b>                     | 6643                                                                                                                                                                                    |
| <b>Gene Symbol</b>                 | SNX2                                                                                                                                                                                    |
| <b>Species</b>                     | Human                                                                                                                                                                                   |
| <b>Immunogen</b>                   | Produced in rabbits immunized with E. coli-derived Human SNX2 fragment                                                                                                                  |
| <b>Notes</b>                       | This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. |
| <b>Product Application Details</b> |                                                                                                                                                                                         |
| <b>Applications</b>                | Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                             |
| <b>Recommended Dilutions</b>       | Immunohistochemistry, Immunocytochemistry/Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin 1:50-1:200                                                                       |




## Images


Immunofluorescence staining of SNX2 in A431 cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with NBP3-21431 (dilution ratio 1:100) at 4C overnight. Then cells were stained with the Alexa Fluor 488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Cytoplasm.



Immunochemical staining of SNX2 in human kidney formalin-fixed paraffin embedded sections with NBP3-21431 at 1:100 dilution



Immunochemical staining of SNX2 in human lymph node formalin-fixed paraffin embedded sections with NBP3-21431 at 1:100 dilution





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-21431-100ul**

---

|                   |                                               |
|-------------------|-----------------------------------------------|
| HAF008            | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7156            | Goat anti-Rabbit IgG (H+L) Secondary Antibody |
| NBP2-24891        | Rabbit IgG Isotype Control                    |
| H00006643-P01-2ug | Recombinant Human SNX2 GST (N-Term) Protein   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-21431](http://www.novusbio.com/reviews/submit/NBP3-21431)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

