

Product Datasheet

C7orf26 Antibody **NBP3-21439-100ul**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-21439

Updated 5/8/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-21439



NBP3-21439-100ul**C7orf26 Antibody**

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.03% Proclin 300
Isotype	IgG
Purity	Antigen and protein A Affinity-purified
Buffer	PBS, pH 7.0


Product Description	
Host	Rabbit
Gene ID	79034
Gene Symbol	C7orf26
Species	Human
Immunogen	Produced in rabbits immunized with E. coli-derived Human C7orf26 fragment
Notes	This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C.

Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Immunohistochemistry, Immunoprecipitation 5-10 uL/mg of lysate, Immunohistochemistry-Paraffin 1:200-1:1000




Images


Immunochemical staining of C7orf26 in human liver formalin-fixed paraffin embedded sections with NBP3-21439 at 1:500 dilution



Immunochemical staining of C7orf26 in human kidney formalin-fixed paraffin embedded sections with NBP3-21439 at 1:500 dilution



C7orf26 was immunoprecipitated using NBP3-21439 at 1:100 dilution, 60 ug of Immunomagnetic beads Protein A/G, HRP Detection Reagent at 1:1000 dilution: Lane A: 0.5 mg RT-4 Whole Cell Lysate. Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 50 kDa. Observed band size: 50 kDa





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-21439-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7156	Goat anti-Rabbit IgG (H+L) Secondary Antibody
NBP2-24891	Rabbit IgG Isotype Control
NBP2-14764PEP	C7orf26 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-21439

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

