

# Product Datasheet

## MVD Antibody (SR1738) NBP3-22438-50ul

Unit Size: 50 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-22438](http://www.novusbio.com/NBP3-22438)

Updated 4/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-22438](http://www.novusbio.com/reviews/destination/NBP3-22438)



**NBP3-22438-50ul**

MVD Antibody (SR1738)


| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 50 ul                                                                                       |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | SR1738                                                                                      |
| Preservative            | 0.02% Sodium Azide                                                                          |
| Isotype                 | IgG                                                                                         |
| Purity                  | Affinity purified                                                                           |
| Buffer                  | PBS, pH 7.4, 150mM NaCl, 50% glycerol.                                                      |
| Target Molecular Weight | 43 kDa                                                                                      |

| Product Description     |                                                                  |
|-------------------------|------------------------------------------------------------------|
| Host                    | Rabbit                                                           |
| Gene Symbol             | MVD                                                              |
| Species                 | Human                                                            |
| Specificity/Sensitivity | Detects endogenous levels of total MVD                           |
| Immunogen               | A synthesized peptide derived from human MVD (Uniprot #: P53602) |

| Product Application Details |                           |
|-----------------------------|---------------------------|
| Applications                | Western Blot              |
| Recommended Dilutions       | Western Blot 1:500-1:2000 |

**Images**

Western Blot: MVD Antibody (SR1738) [NBP3-22438] - Western blot analysis using NBP3-22438 in K562 cell lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-22438-50ul**

---

|                 |                                                     |
|-----------------|-----------------------------------------------------|
| HAF008          | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160          | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891      | Rabbit IgG Isotype Control                          |
| NBP1-78853-50ug | Recombinant Human MVD His Protein                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-22438](http://www.novusbio.com/reviews/submit/NBP3-22438)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

