

Product Datasheet

Integrin beta 4/CD104 [p Tyr1526] Antibody NBP3-23401

Unit Size: 0.1 ml

Store at -20C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-23401

Updated 7/12/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-23401



NBP3-23401**Integrin beta 4/CD104 [p Tyr1526] Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1 mg/ml BSA and 50% glycerol
Product Description	
Description	Stable for at least 1 year at -20°C.
Host	Rabbit
Gene ID	3691
Gene Symbol	ITGB4
Species	Human, Mouse, Rat
Specificity/Sensitivity	This antibody was cross-adsorbed to phospho-tyrosine coupled to agarose and to unphosphorylated Integrin beta 4/CD104 (Tyr-1526) peptide before affinity purification using phospho-Integrin beta 4/CD104 (Tyr-1526) peptide (without carrier). The purified antibody detects a 200kDa protein corresponding to the molecular mass of Integrin beta 4/CD104 on SDS-PAGE immunoblots of A431 and src-transformed Hct116 cells treated with pervanadate, but not in control cells.
Immunogen	Phospho-Integrin beta 4/CD104 (Tyr-1526) synthetic peptide (coupled to KLH) corresponding to amino acid residues around tyrosine 1526 of human Integrin beta 4/CD104. This peptide sequence is found in all three Integrin beta 4/CD104 isoforms and is highly conserved in rat and mouse Integrin beta 4/CD104. [UniProt# P16144]
Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1:1000, Immunocytochemistry/Immunofluorescence 1:100





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-23401

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7156	Goat anti-Rabbit IgG (H+L) Secondary Antibody
NBP2-24891	Rabbit IgG Isotype Control
NBP2-61461-10ug	Recombinant Human Integrin beta 4/CD104 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-23401

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

