

# Product Datasheet

## Androgen Receptor V7 Antibody (DHTR.V7/9394R) [Janelia Fluor® 635] NBP3-24295JF635

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-24295JF635](http://www.novusbio.com/NBP3-24295JF635)

Updated 8/19/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-24295JF635](http://www.novusbio.com/reviews/destination/NBP3-24295JF635)



**NBP3-24295JF635**

Androgen Receptor V7 Antibody (DHTR.V7/9394R) [Janelia Fluor® 635]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Clone                       | DHTR.V7/9394R                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Isotype                     | IgG Kappa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Conjugate                   | Janelia Fluor 635                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Purity                      | Protein A purified                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Gene ID                     | 367                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Gene Symbol                 | AR                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Specificity/Sensitivity     | AR V7 belongs to the nuclear hormone receptor family. It contains the nuclear receptor DNA binding domain. Steroid hormone receptors are ligand-activated transcription factors that regulate the expression of eukaryotic genes affecting cell proliferation and differentiation in target tissues. The activity of transcription factors is regulated by binding co-activators and co-repressors. This target is a splice variant of the androgen receptor that lacks the C-terminal androgen-binding site. AR-V7 plays an important role in the occurrence, development and drug resistance of prostate cancer. |
| Immunogen                   | Recombinant fragment (around aa600-920) of human Androgen Receptor-protein (Androgen Receptor-V7 specific) (exact sequence is proprietary)                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Applications                | Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions       | Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-24295JF635](http://www.novusbio.com/reviews/submit/NBP3-24295JF635)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

