

# Product Datasheet

## DLG7/HURP Recombinant Protein Antigen NBP3-24790PEP

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-24790PEP](http://www.novusbio.com/NBP3-24790PEP)

Updated 2/3/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-24790PEP](http://www.novusbio.com/reviews/destination/NBP3-24790PEP)



**NBP3-24790PEP****DLG7/HURP Recombinant Protein Antigen****Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| <b>Preservative</b>  | No Preservative                                                                             |
| <b>Purity</b>        | >80% by SDS-PAGE and Coomassie blue staining                                                |
| <b>Buffer</b>        | PBS and 1M Urea, pH 7.4.                                                                    |

**Product Description**

|                    |                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | <p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human DLG7/HURP</p> <p><b>Source:</b> <i>E.coli</i></p> <p><b>Amino Acid Sequence:</b><br/>QTRLLMKERFKQFEGLVDDCEYKRGIKETTCTDLDGFWDMSFQIEDVIHKFN<br/>NLIKLEESGWQVNNNMNHNMNKNVFRKKVVSGIASKPKQDDAGRI</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p> |
| <b>Gene ID</b>     | 9787                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Gene Symbol</b> | DLGAP5                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Species</b>     | Human                                                                                                                                                                                                                                                                                                                                                                                                    |

**Product Application Details**

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Antibody Competition                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Recommended Dilutions</b> | Antibody Competition 10-100 molar excess                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Application Notes</b>     | <p>This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP3-24790. It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml. For current lot information, including availability, please contact our technical support team click <a href="mailto:nb-technical@bio-techne.com">nb-technical@bio-techne.com</a>.</p> |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-24790PEP](http://www.novusbio.com/reviews/submit/NBP3-24790PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

