

# Product Datasheet

## KCNA7 Antibody - BSA Free

### NBP3-24934

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-24934](http://www.novusbio.com/NBP3-24934)

Updated 2/25/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-24934](http://www.novusbio.com/reviews/destination/NBP3-24934)



**NBP3-24934**

KCNA7 Antibody - BSA Free

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                                               |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Affinity purified                                                                                                |
| <b>Buffer</b>        | PBS, pH 7.2, 40% glycerol                                                                                        |

**Product Description**


|                         |                                                                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>             | Rabbit                                                                                                                                                                  |
| <b>Gene ID</b>          | 3743                                                                                                                                                                    |
| <b>Gene Symbol</b>      | KCNA7                                                                                                                                                                   |
| <b>Species</b>          | Human                                                                                                                                                                   |
| <b>Reactivity Notes</b> | Mouse, Rat (88%).                                                                                                                                                       |
| <b>Immunogen</b>        | This antibody has been engineered to specifically recognize the recombinant protein KCNA7 using the following amino acid sequence:<br>TEGEEAGMFSHVDMQPCGPLEGKANGGLVDGEV |

**Product Application Details**

|                              |                                                                                                                                    |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Immunocytochemistry/ Immunofluorescence                                                                                            |
| <b>Recommended Dilutions</b> | Immunocytochemistry/ Immunofluorescence 0.25-2 µg/ml                                                                               |
| <b>Application Notes</b>     | For ICC/IF, we recommend using a combination of PFA and Triton X-100. This will give you the optimal results for your experiments. |

**Images**

Staining of human cell line SiHa shows localization to vesicles.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-24934**

---

|            |                                                     |
|------------|-----------------------------------------------------|
| HAF008     | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160     | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control                          |
| 355-BM-010 | BMP-2 [Unconjugated]                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-24934](http://www.novusbio.com/reviews/submit/NBP3-24934)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

