

# Product Datasheet

## African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free NBP3-25367

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-25367](http://www.novusbio.com/NBP3-25367)

Updated 7/31/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-25367](http://www.novusbio.com/reviews/destination/NBP3-25367)



**NBP3-25367**

African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1219
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	22220355
Gene Symbol	E183L
Species	Virus
Reactivity Notes	African Swine Fever Virus
Immunogen	Recombinant protein encompassing a sequence within the center region of African Swine Fever Virus p54. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA, Lateral Flow Assay
Recommended Dilutions	Western Blot 1:1000-1:10000, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA, Lateral Flow Assay

**Images**


ELISA: African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free [NBP3-25367] - Sandwich ELISA detection of recombinant ASFV p54 using antibodies as below. Capture: ASFV p54 antibody [HL1218] (NBP3-25366) (5 ug/mL) Detection: HRP-conjugated ASFV p54 antibody [HL1219] (NBP3-25367) (0.1 ug/mL) Please notice that NBP3-25367 needs to be conjugated to HRP to function as the detection antibody when paired with NBP3-25366. Please contact us for custom HRP-conjugated antibody.




Immunocytochemistry/Immunofluorescence: African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free [NBP3-25367] - ASFV p54 antibody [HL1219] detects ASFV p54 protein by immunofluorescent analysis. Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: ASFV p54 stained by ASFV p54 antibody [HL1219] (NBP3-25367) diluted at 1:500. Blue: Fluoroshield with DAPI . Scale bar= 10um.




Western Blot: African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free [NBP3-25367] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with ASFV p54 antibody [HL1219] (NBP3-25367) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




ELISA: African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free [NBP3-25367] - Indirect ELISA analysis performed by coating plate with recombinant recombinant HEK293 expressed, full-length ASFV p54 (40-0.63 nM). Coated protein was probed with ASFV p54 antibody [HL1219] (NBP3-25367) (1 ug/mL). Goat anti-rabbit IgG antibody (HRP) (1:10000) was used to detect bound primary antibody.



Lateral Flow Assay: African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free [NBP3-25367] - Detection of HEK293 expressed, ASFV p54 (ECD) protein by lateral flow assay using the monoclonal antibody pair. Capture: ASFV p54 antibody (NBP3-25367 [HL1219]) Detection: ASFV p54 antibody (NBP3-25370 [HL1289]) Samples (20 ng) : 1. ASFV p54 (ECD) protein 2. Lysis buffer



Lateral Flow Assay: African Swine Fever Virus p54 Antibody (HL1219) - Azide and BSA Free [NBP3-25367] - Detection of HEK293 expressed, ASFV p54 (ECD) protein by lateral flow assay using the monoclonal antibody pair. Capture: ASFV p54 antibody (NBP3-25367 [HL1219])  
Detection: ASFV p54 antibody (NBP3-25366 [HL1218]) Samples (20 ng)  
: 1. ASFV p54 (ECD) protein 2. Lysis buffer





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-25367**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-25367](http://www.novusbio.com/reviews/submit/NBP3-25367)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

