

Product Datasheet

ATP6V1H Antibody (HL1316) - Azide and BSA Free NBP3-25372

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25372

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25372



NBP3-25372

ATP6V1H Antibody (HL1316) - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1316
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description


Host	Rabbit
Gene ID	51606
Gene Symbol	ATP6V1H
Species	Human, Mouse, Drosophila
Immunogen	Recombinant protein encompassing a sequence within the center region of human ATP6V1H. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry, Immunohistochemistry-Paraffin

Images


Western Blot: ATP6V1H Antibody (HL1316) - Azide and BSA Free [NBP3-25372] - Drosophila tissue extract (50 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with ATP6V1H antibody [HL1316] (NBP3-25372) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunohistochemistry-Paraffin: ATP6V1H Antibody (HL1316) - Azide and BSA Free [NBP3-25372] - ATP6V1H antibody [HL1316] detects ATP6V1H protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse kidney. ATP6V1H stained by ATP6V1H antibody [HL1316] (NBP3-25372) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: ATP6V1H Antibody (HL1316) - Azide and BSA Free [NBP3-25372] - Various whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with ATP6V1H antibody [HL1316] (NBP3-25372) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunohistochemistry-Paraffin: ATP6V1H Antibody (HL1316) - Azide and BSA Free [NBP3-25372] - ATP6V1H antibody [HL1316] detects ATP6V1H protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse heart. ATP6V1H stained by ATP6V1H antibody [HL1316] (NBP3-25372) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: ATP6V1H Antibody (HL1316) - Azide and BSA Free [NBP3-25372] - ATP6V1H antibody [HL1316] detects ATP6V1H protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded mouse muscle. ATP6V1H stained by ATP6V1H antibody [HL1316] (NBP3-25372) diluted at 1:200. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: ATP6V1H Antibody (HL1316) - Azide and BSA Free [NBP3-25372] - ATP6V1H antibody [HL1316] detects ATP6V1H protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissues. ATP6V1H stained by ATP6V1H antibody [HL1316] (NBP3-25372) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25372

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00051606-P01-10ug	Recombinant Human ATP6V1H GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25372

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

