

Product Datasheet

CXCL13/BLC/BCA-1 Antibody (HL2172) - Azide and BSA Free **NBP3-25379**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25379

Updated 7/31/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25379



NBP3-25379

CXCL13/BLC/BCA-1 Antibody (HL2172) - Azide and BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2172
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS


Product Description	
Host	Rabbit
Gene ID	10563
Gene Symbol	CXCL13
Species	Human
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human CXCL13/BLC/BCA-1. The exact sequence is proprietary.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry, Immunohistochemistry-Paraffin




Images


Immunohistochemistry-Paraffin: CXCL13/BLC/BCA-1 Antibody (HL2172) - Azide and BSA Free [NBP3-25379] - BCA1 antibody [HL2172] detects BCA1 protein by immunohistochemical analysis. Sample: Paraffin-embedded human Clear Cell Renal Cell Carcinoma (ccRCC). BCA1 stained by BCA1 antibody [HL2172] (NBP3-25379) diluted at 1:100. Antigen Retrieval: Tris-EDTA buffer, pH 8.0, 15 min



Immunohistochemistry-Paraffin: CXCL13/BLC/BCA-1 Antibody (HL2172) - Azide and BSA Free [NBP3-25379] - BCA1 antibody [HL2172] detects secreted BCA1 protein by immunohistochemical analysis. Sample: Paraffin-embedded human spleen. BCA1 stained by BCA1 antibody [HL2172] (NBP3-25379) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: CXCL13/BLC/BCA-1 Antibody (HL2172) - Azide and BSA Free [NBP3-25379] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with BCA1 antibody [HL2172] (NBP3-25379) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25379

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-34885-5ug	Recombinant Human CXCL13/BLC/BCA-1 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25379

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

