

Product Datasheet

CCR4 Antibody (HL2481) - Azide and BSA Free NBP3-25402

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25402

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25402



NBP3-25402

CCR4 Antibody (HL2481) - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL2481
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS
Target Molecular Weight	70 kDa

Product Description


Host	Rabbit
Gene ID	1233
Gene Symbol	CCR4
Species	Human, Mouse, Rat
Immunogen	Synthetic peptide of human CCR4.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry Assay dependent, Immunocytochemistry/ Immunofluorescence Assay dependent, Immunohistochemistry-Paraffin Assay dependent

Images


Western Blot: CCR4 Antibody (HL2481) - Azide and BSA Free [NBP3-25402] - Non-transfected (-) and transfected (+) boiled and unboiled 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with CCR4 antibody [HL2481] (NBP3-25402) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunohistochemistry-Paraffin: CCR4 Antibody (HL2481) - Azide and BSA Free [NBP3-25402] - CCR4 antibody [HL2481] detects CCR4 protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded rat thymus gland. CCR4 stained by CCR4 antibody [HL2481] (NBP3-25402) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Western Blot: CCR4 Antibody (HL2481) - Azide and BSA Free [NBP3-25402] - Mouse tissue extract (50 ug) was separated by 7.5% SDS-PAGE, and the membrane was blotted with CCR4 antibody [HL2481] (NBP3-25402) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Western Blot: CCR4 Antibody (HL2481) - Azide and BSA Free [NBP3-25402] - Non-transfected (-) and transfected (+) boiled and unboiled SW480 whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with CCR4 antibody [HL2481] (NBP3-25402) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Western Blot: CCR4 Antibody (HL2481) - Azide and BSA Free [NBP3-25402] - Unboiled SW480 whole cell and membrane extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with CCR4 antibody [HL2481] (NBP3-25402) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced. (WCE: whole cell extract; ME: membrane extract)



Immunohistochemistry-Paraffin: CCR4 Antibody (HL2481) - Azide and BSA Free [NBP3-25402] - CCR4 antibody [HL2481] detects CCR4 protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded mouse thymus gland. CCR4 stained by CCR4 antibody [HL2481] (NBP3-25402) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



ICC/IF: CCR4 Antibody (HL2481) -SW480 cells were fixed in ice-cold MeOH for 5 min. Green: CCR4 stained by CCR4 antibody [HL2481] diluted at 1:200. Blue: Fluoroshield with DAPI.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25402

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB100-56336PEP	CCR4 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25402

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

