

Product Datasheet

FCER1G Antibody (HL1418) - Azide and BSA Free NBP3-25463

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25463

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25463



NBP3-25463

FCER1G Antibody (HL1418) - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1418
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description


Host	Rabbit
Gene ID	2207
Gene Symbol	FCER1G
Species	Human, Mouse, Rat, Canine, Feline
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human FCER1G. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin

Images


Immunohistochemistry-Paraffin: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - FCER1G antibody [HL1418] detects FCER1G protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded dog cerebellum. FCER1G stained by FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Western Blot: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - Whole cell extract (30 ug) was separated by 15% SDS-PAGE, and the membrane was blotted with FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Western Blot: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunohistochemistry-Paraffin: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - FCER1G antibody [HL1418] detects FCER1G protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded cat cerebellum. FCER1G stained by FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Immunocytochemistry/Immunofluorescence: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - FCER1G antibody [HL1418] detects FCER1G protein at cell membrane by immunofluorescent analysis. Sample: THP-1 cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: FCER1G stained by FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:500. Blue: Fluoroshield with DAPI .




Immunohistochemistry-Paraffin: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - FCER1G antibody [HL1418] detects FCER1G protein by immunohistochemical analysis. Sample: Paraffin-embedded rat tissues. FCER1G stained by FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - FCER1G antibody [HL1418] detects FCER1G protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded mouse spleen. FCER1G stained by FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: FCER1G Antibody (HL1418) - Azide and BSA Free [NBP3-25463] - FCER1G antibody [HL1418] detects FCER1G protein at cell membrane by immunohistochemical analysis. Sample: Paraffin-embedded mouse spleen. FCER1G stained by FCER1G antibody [HL1418] (NBP3-25463) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25463

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00002207-P01-10ug	Recombinant Human FCER1G GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25463

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

