

# Product Datasheet

## Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free NBP3-25524

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-25524](http://www.novusbio.com/NBP3-25524)

Updated 7/31/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-25524](http://www.novusbio.com/reviews/destination/NBP3-25524)



**NBP3-25524**

Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1273
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Host	Rabbit
Species	Virus
Immunogen	Recombinant fragment corresponding to C-terminal region of Influenza A virus M1 (A/Puerto Rico/8/1934(H1N1)). The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin

**Images**


Immunocytochemistry/Immunofluorescence: Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free [NBP3-25524] - Immunofluorescent analysis of Influenza virus infected cells using Influenza A virus M1 (matrix protein) antibody [HL1273] (NBP3-25524). Sample: Influenza A and B virus infected cells slide. Green: Influenza A virus M1 (matrix protein) antibody [HL1273] (NBP3-25524) diluted at 1:100. Blue: Fluoroshield with DAPI .




**Western Blot: Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free [NBP3-25524]** - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with Influenza A virus M1 (matrix protein) antibody [HL1273] (NBP3-25524) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




**Western Blot: Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free [NBP3-25524]** - Non-infected (-) and infected (+) H1299 whole cell extracts (15 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with Influenza A virus M1 (matrix protein) antibody [HL1273] (NBP3-25524) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



**Immunohistochemistry-Paraffin: Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free [NBP3-25524]** - Influenza A virus M1 (matrix protein) antibody [HL1273] detects Influenza A virus M1 (matrix protein) protein at cytoplasm by immunohistochemical analysis. Sample: Paraffin-embedded Influenza A virus M1 transfected 293T. Green: Influenza A virus M1 (matrix protein) stained by Influenza A virus M1 (matrix protein) antibody [HL1273] (NBP3-25524) diluted at 1:1000. Blue: Fluoroshield with DAPI. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



**Immunocytochemistry/Immunofluorescence: Influenza A H1N1 M1 Antibody (HL1273) - Azide and BSA Free [NBP3-25524]** - Influenza A virus M1 (matrix protein) antibody [HL1273] detects Influenza A virus M1 (matrix protein) protein by immunofluorescent analysis. Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: Influenza A virus M1 (matrix protein) stained by Influenza A virus M1 (matrix protein) antibody [HL1273] (NBP3-25524) diluted at 1:500. Blue: Fluoroshield with DAPI.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-25524**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-25524](http://www.novusbio.com/reviews/submit/NBP3-25524)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

