

Product Datasheet

NFIX Antibody (HL1269) - Azide and BSA Free NBP3-25576

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-25576

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-25576



NBP3-25576

NFIIX Antibody (HL1269) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1269
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	4784
Gene Symbol	NFIIX
Species	Mouse, Rat
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Xenopus (87%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human NFIIX. The exact sequence is proprietary.


Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Immunohistochemistry-Paraffin: NFIIX Antibody (HL1269) - Azide and BSA Free [NBP3-25576] - NFIIX antibody [HL1269] detects NFIIX protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissue. NFIIX stained by NFIIX antibody [HL1269] (NBP3-25576) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Immunocytochemistry/Immunofluorescence: NFIX Antibody (HL1269) - Azide and BSA Free [NBP3-25576] - NFIX antibody [HL1269] detects NFIX protein at nucleus by immunofluorescent analysis. Sample: DIV10 rat E18 primary hippocampal neuron and glia cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: NFIX stained by NFIX antibody [HL1269] (NBP3-25576) diluted at 1:250. Red: Tau, a cytoskeleton marker, stained by Tau antibody [GT287] diluted at 1:500.




Immunohistochemistry-Paraffin: NFIX Antibody (HL1269) - Azide and BSA Free [NBP3-25576] - NFIX antibody [HL1269] detects NFIX protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse brain. NFIX stained by NFIX antibody [HL1269] (NBP3-25576) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: NFIX Antibody (HL1269) - Azide and BSA Free [NBP3-25576] - NFIX antibody [HL1269] detects NFIX protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded rat brain. NFIX stained by NFIX antibody [HL1269] (NBP3-25576) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: NFIX Antibody (HL1269) - Azide and BSA Free [NBP3-25576] - NFIX antibody [HL1269] detects NFIX protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse brain. NFIX stained by NFIX antibody [HL1269] (NBP3-25576) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-25576

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-58904PEP	NFIX Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-25576

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

