

Product Datasheet

Somatostatin R2/SSTR2 Antibody (7A9) - BSA Free NBP3-26301-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-26301

Updated 7/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-26301



NBP3-26301-100ul

Somatostatin R2/SSTR2 Antibody (7A9) - BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20 to -70C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	7A9
Preservative	0.03% Proclin 300
Isotype	IgG1
Purity	Affinity purified
Buffer	50% Glycerol, 0.01M PBS, PH 7.4

Product Description

Host	Human
Gene ID	6752
Gene Symbol	SSTR2
Species	Human
Immunogen	Recombinant Human Somatostatin R2/SSTR2 protein [UniProt P30874]

Product Application Details

Applications	ELISA, Flow Cytometry
Recommended Dilutions	Flow Cytometry, ELISA




Images


ELISA: Somatostatin R2/SSTR2 Antibody (7A9) [NBP3-26301] - The Binding Activity of Human Somatostatin R2/SSTR2 with Somatostatin R2/SSTR2 Antibody.

Activity: Measured by its binding ability in a functional ELISA.

Immobilized Human Somatostatin R2/SSTR2 at 10 ug/mL can bind Somatostatin R2/SSTR2 Antibody, the EC₅₀ is 58.13-81.28 ng/mL.



Flow Cytometry: Somatostatin R2/SSTR2 Antibody (7A9) [NBP3-26301]
 - Untransfected CT26 cells (green line) and transfected Human SSTR2 CT26 stable cells (red line) were stained with anti-SSTR2 recombinant antibody (2ug/1*10⁶ cells), washed and then followed by APC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-26301-100ul

G-102-C	Goat anti-Human IgG Secondary Antibody [Unconjugated]
NB7446	Goat anti-Human IgG Fc Secondary Antibody
DDXCH01P-100	Human IgG1 Isotype Control
NB300-157PEP	Somatostatin R2/SSTR2 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-26301

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

