

# Product Datasheet

## Uteroglobin/SCGB1A1 Antibody (21B2) - BSA Free NBP3-26335-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-26335](http://www.novusbio.com/NBP3-26335)

Updated 7/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-26335](http://www.novusbio.com/reviews/destination/NBP3-26335)



**NBP3-26335-100ul**

Uteroglobin/SCGB1A1 Antibody (21B2) - BSA Free

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at -20 to -70C. Avoid freeze-thaw cycles.                                             |
| <b>Clonality</b>     | Monoclonal                                                                                  |
| <b>Clone</b>         | 21B2                                                                                        |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Affinity purified                                                                           |
| <b>Buffer</b>        | PBS, pH 7.4, 150mM NaCl, and 50% glycerol                                                   |

**Product Description**


|                    |                                                                               |
|--------------------|-------------------------------------------------------------------------------|
| <b>Host</b>        | Rabbit                                                                        |
| <b>Gene ID</b>     | 7356                                                                          |
| <b>Gene Symbol</b> | SCGB1A1                                                                       |
| <b>Species</b>     | Human                                                                         |
| <b>Immunogen</b>   | A synthesized peptide derived from Human Uteroglobin/SCGB1A1 [UniProt P11684] |

**Product Application Details**

|                              |                                  |
|------------------------------|----------------------------------|
| <b>Applications</b>          | ELISA, Flow Cytometry            |
| <b>Recommended Dilutions</b> | Flow Cytometry 1:50-1:200, ELISA |

**Images**

Flow Cytometry: Uteroglobin/SCGB1A1 Antibody (21B2) [NBP3-26335] - Overlay Peak curve showing A549 cells stained with Uteroglobin/SCGB1A1 Antibody (21B2) (red line) at 1:50. The cells were fixed in 4% formaldehyde and permeated by 0.2% Triton X-100. Then 10% normal goat serum to block non-specific protein-protein interactions followed by the antibody (1ug/1\*10<sup>6</sup> cells) for 45min at 4C. The secondary antibody used was FITC-conjugated Goat Anti-rabbit IgG (H +L) at 1:200 dilution for 35min at 4C. Control antibody (green line) was rabbit IgG (1ug/1\*10<sup>6</sup> cells) used under the same conditions. Acquisition of >10,000 events was performed.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-26335-100ul**

---

|                 |                                                     |
|-----------------|-----------------------------------------------------|
| HAF008          | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160          | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891      | Rabbit IgG Isotype Control                          |
| NBP2-35126-10ug | Recombinant Human Uteroglobulin/SCGB1A1 Protein     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-26335](http://www.novusbio.com/reviews/submit/NBP3-26335)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

