

# Product Datasheet

## NSDHL Antibody (20F2) - BSA Free NBP3-26457-100ul

Unit Size: 100 ul

Store at -20 to -70C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-26457](http://www.novusbio.com/NBP3-26457)

Updated 7/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-26457](http://www.novusbio.com/reviews/destination/NBP3-26457)



**NBP3-26457-100ul**


NSDHL Antibody (20F2) - BSA Free

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at -20 to -70C. Avoid freeze-thaw cycles.                                             |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | 20F2                                                                                        |
| Preservative                | 0.02% Sodium Azide                                                                          |
| Isotype                     | IgG                                                                                         |
| Purity                      | Affinity purified                                                                           |
| Buffer                      | PBS, pH 7.4, 150mM NaCl, and 50% glycerol                                                   |
| Product Description         |                                                                                             |
| Host                        | Rabbit                                                                                      |
| Gene ID                     | 50814                                                                                       |
| Gene Symbol                 | NSDHL                                                                                       |
| Species                     | Human                                                                                       |
| Immunogen                   | A synthesized peptide derived from Human NSDHL [UniProt Q15738]                             |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, ELISA, Immunohistochemistry                                                   |
| Recommended Dilutions       | Western Blot 1:500-1:2000, ELISA, Immunohistochemistry 1:50-1:200                           |



## Images

**Immunohistochemistry:** NSDHL Antibody (20F2) [NBP3-26457] - Image of NSDHL Antibody (20F2) diluted at 1:100 and staining in paraffin-embedded human testis tissue performed. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4C overnight. The primary is detected by a Goat anti-rabbit polymer IgG labeled by HRP and visualized using 0.05% DAB.




**Western Blot:** NSDHL Antibody (20F2) [NBP3-26457] - Positive Western Blot detected in: HepG2 whole cell lysate, Hela whole cell lysate, PC-3 whole cell lysate, HEK293 whole cell lysate.

All lanes: NSDHL Antibody at 1: 1000

Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution.

Predicted band size: 42 kDa

Observed band size: 42 kDa





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-26457-100ul**

---

|            |                                                     |
|------------|-----------------------------------------------------|
| HAF008     | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160     | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control                          |
| NBP2-23258 | Recombinant Human NSDHL His Protein                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-26457](http://www.novusbio.com/reviews/submit/NBP3-26457)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

