

Product Datasheet

SMCX Antibody (PCRP-KDM5C-1A11) - Azide and BSA Free NBP3-27035

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27035

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27035



NBP3-27035

SMCX Antibody (PCRP-KDM5C-1A11) - Azide and BSA Free

| Product Information | |
|---------------------|---|
| Unit Size | 100 ug |
| Concentration | 1 mg/ml |
| Storage | Store at -20 to -80C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | PCRP-KDM5C-1A11 |
| Preservative | No Preservative |
| Isotype | IgG2c |
| Purity | Protein A or G purified |
| Buffer | 10mM PBS |

| Product Description | |
|---------------------|---|
| Description | Antibody without azide - store at -20 to -80C. Non-hazardous. No MSDS required. |
| Host | Mouse |
| Gene ID | 8242 |
| Gene Symbol | KDM5C |
| Species | Human |
| Immunogen | Recombinant full-length human SMCX protein |

| Product Application Details | |
|-----------------------------|---|
| Applications | Flow Cytometry |
| Recommended Dilutions | Flow Cytometry 1-2ug/million cells |
| Application Notes | Optimal dilution for a specific application should be determined. |



Images

Analysis of Protein Array containing more than 19,000 full-length human proteins using KDM5C Mouse Monoclonal (PCRP-KDM5C-1A11).



Flow Cytometric Analysis of PFA-fixed HeLa cells. SMCX Antibody (PCRP-KDM5C-1A11) - Azide and BSA Free followed by goat anti-mouse IgG-CF488 (blue); unstained cells (red).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27035

| | |
|---------------|---|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB720-B | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP2-55541PEP | SMCX Recombinant Protein Antigen |
| E-305-025 | Ubiquitin-activating Enzyme/UBE1 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27035

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

