

Product Datasheet

TLX3 Antibody - BSA Free NBP3-27888

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27888

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27888



NBP3-27888

TLX3 Antibody - BSA Free

Product Information

| | |
|----------------------|--|
| Unit Size | 0.1 mg |
| Concentration | 1 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Purity | Peptide affinity purified |
| Buffer | PBS |

Product Description

| | |
|--------------------------------|---|
| Description | Novus Biologicals Rabbit TLX3 Antibody - BSA Free (NBP3-27888) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Rabbit |
| Gene ID | 30012 |
| Gene Symbol | TLX3 |
| Species | Human, Mouse, Rat |
| Specificity/Sensitivity | TLX3 antibody is predicted to not cross-react with other TLX proteins. |
| Immunogen | TLX3 antibody was raised against a 16 amino acid peptide near the amino terminus of human TLX3. The immunogen is located within the first 50 amino acids of TLX3. |


Product Application Details

| | |
|------------------------------|---|
| Applications | Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence |
| Recommended Dilutions | Western Blot 0.5 - 1 ug/mL, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin |




Images


Western Blot: TLX3 Antibody [NBP3-27888] - Western blot analysis of TLX3 in human spleen tissue lysate with TLX3 antibody at (A) 0.5 and (B) 1 u/mL.



Immunohistochemistry: TLX3 Antibody [NBP3-27888] - Immunohistochemistry of TLX3 in human spleen tissue with TLX3 antibody at 2.5 u/ml.



Immunocytochemistry/ Immunofluorescence: TLX3 Antibody [NBP3-27888] - Immunofluorescence of TLX3 in human spleen tissue with TLX3 antibody at 20 u/ml.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27888

| | |
|--------------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| H00030012-P01-10ug | Recombinant Human TLX3 GST (N-Term) Protein |
| NBL1-16958 | TLX3 Overexpression Lysate |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27888

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

