

Product Datasheet

PEX3 Antibody - BSA Free NBP3-27906

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-27906

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-27906



NBP3-27906

PEX3 Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Purity	Peptide affinity purified
Buffer	PBS

Product Description


Description	Novus Biologicals Rabbit PEX3 Antibody - BSA Free (NBP3-27906) is a polyclonal antibody validated for use in WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8504
Gene Symbol	PEX3
Species	Human, Mouse, Rat
Immunogen	PEX3 antibody was raised against a peptide corresponding to 14 amino acids near the center of human PEX3. The immunogen is located within the center of PEX3.

Product Application Details


Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 2 ug/mL, ELISA, Immunocytochemistry/ Immunofluorescence 20 ug/mL

Images


Immunocytochemistry/ Immunofluorescence: PEX3 Antibody [NBP3-27906] - Immunofluorescent analysis of 4% paraformaldehyde-fixed human spleen labeling PEX3 at 20 u/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).




Western Blot: PEX3 Antibody [NBP3-27906] - Loading: 15 u of lysate
 Antibodies: PEX3, 2 u/mL, 1h incubation at RT in 5% NFDM/TBST.
 Secondary: Goat Anti-Rabbit IgG HRP conjugate at 1:10000 dilution.




Immunocytochemistry/ Immunofluorescence: PEX3 Antibody [NBP3-27906] - Immunofluorescent analysis of 4% paraformaldehyde-fixed mouse thymus labeling PEX3 at 20 u/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).




Immunocytochemistry/ Immunofluorescence: PEX3 Antibody [NBP3-27906] - Immunofluorescent analysis of 4% paraformaldehyde-fixed rat testis labeling PEX3 at 20 u/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).




Immunocytochemistry/ Immunofluorescence: PEX3 Antibody [NBP3-27906] - Immunofluorescent analysis of 4% paraformaldehyde-fixed HeLa cells labeling PEX3 at 20 u/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).



Western Blot: PEX3 Antibody [NBP3-27906] - Loading: 15 u of lysates per lane. Antibodies: PEX3, 2 u/mL, 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Western Blot: PEX3 Antibody [NBP3-27906] - Loading: 15 u of lysates per lane. Antibodies: PEX3, 2 u/mL, 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-27906

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-38838PEP	PEX3 Recombinant Protein Antigen
8090-ZN-025	PHEX

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-27906

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

