

# Product Datasheet

## RSU1 Antibody - BSA Free NBP3-28946

Unit Size: 100 ug

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-28946](http://www.novusbio.com/NBP3-28946)

Updated 3/18/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-28946](http://www.novusbio.com/reviews/destination/NBP3-28946)



**NBP3-28946**


RSU1 Antibody - BSA Free

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ug                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C. Do not freeze.                                                                 |
| Clonality                   | Polyclonal                                                                                  |
| Preservative                | 0.09% Sodium Azide                                                                          |
| Isotype                     | IgG                                                                                         |
| Purity                      | Affinity purified                                                                           |
| Buffer                      | Tris-citrate/phosphate pH 7 -8.                                                             |
| Product Description         |                                                                                             |
| Host                        | Rabbit                                                                                      |
| Gene ID                     | 6251                                                                                        |
| Gene Symbol                 | RSU1                                                                                        |
| Species                     | Human, Mouse                                                                                |
| Immunogen                   | Between 227 and 277                                                                         |
| Assay Type                  | [Assay Type]                                                                                |
| Sample Volume               | [Sample Volume]                                                                             |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, Immunoprecipitation                                                           |
| Recommended Dilutions       | Western Blot, Immunoprecipitation                                                           |




## Images

Immunoprecipitation:[NBP3-28946]-Detection of human RSU1 by western blot of immunoprecipitates. Samples: Whole cell lysate from 293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-RSU1 antibody used for IP at 6  $\mu$ g per reaction. RSU1 was also immunoprecipitated by rabbit anti-RSU1 antibody. For blotting immunoprecipitated RSU1, was used at 0.4  $\mu$ g/ml. Detection: Chemiluminescence with an exposure time of 10 seconds.



Western Blot:[NBP3-28946]-RSU1 Antibody





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-28946**

|                  |                                                     |
|------------------|-----------------------------------------------------|
| HAF008           | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160           | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891       | Rabbit IgG Isotype Control                          |
| NBP2-51675-0.1mg | Recombinant Human RSU1 His Protein                  |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-28946](http://www.novusbio.com/reviews/submit/NBP3-28946)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

