

# Product Datasheet

## CBLL1 Antibody - BSA Free

### NBP3-30247

Unit Size: 100 ug

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-30247](http://www.novusbio.com/NBP3-30247)

Updated 3/18/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-30247](http://www.novusbio.com/reviews/destination/NBP3-30247)



**NBP3-30247**

## CBLL1 Antibody - BSA Free

## Product Information

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ug                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C. Do not freeze.                                                                 |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.09% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Affinity purified                                                                           |
| <b>Buffer</b>        | Tris-citrate/phosphate buffer, pH 7 to 8                                                    |

## Product Description


|               |                     |
|---------------|---------------------|
| Host          | Rabbit              |
| Gene ID       | 79872               |
| Gene Symbol   | CBLL1               |
| Species       | Human, Mouse        |
| Immunogen     | between 441 and 491 |
| Assay Type    | [Assay Type]        |
| Sample Volume | [Sample Volume]     |

## Product Application Details

|                              |                                   |
|------------------------------|-----------------------------------|
| <b>Applications</b>          | Western Blot, Immunoprecipitation |
| <b>Recommended Dilutions</b> | Western Blot, Immunoprecipitation |

## Images

**Western Blot:**[NBP3-30247]-Detection of human and mouse Hakai by western blot and immunoprecipitation. Samples: Whole cell lysate from HeLa , HEK293T and mouse NIH 3T3 cells. Antibodies: Affinity purified rabbit anti-Hakai antibody used for WB at 0. 1 µg/ml and 1 µg/ml and used for IP at 3 µg/mg lysate. Hakai was also immunoprecipitated by rabbit anti-Hakai antibody , which recognizes an upstream epitope. Detection: Chemiluminescence with exposure times of 30 seconds and 10 seconds.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-30247**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-83589PEP | CBLL1 Recombinant Protein Antigen                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-30247](http://www.novusbio.com/reviews/submit/NBP3-30247)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

