

Product Datasheet

EXOC2 Antibody (JE65-98)

NBP3-32322

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-32322

Updated 8/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-32322



NBP3-32322

EXOC2 Antibody (JE65-98)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	JE65-98
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	1*TBS (pH7.4), 0.05% BSA and 40% Glycerol
Target Molecular Weight	104 kDa
Product Description	
Host	Rabbit
Gene ID	55770
Gene Symbol	EXOC2
Species	Human, Mouse, Rat
Immunogen	Recombinant protein within Human EXOC2 aa 1-100 / 924. (Uniprot: Q96KP1)
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry, Immunohistochemistry-Paraffin 1:200



Images

Western Blot: EXOC2 Antibody (JE65-98) [NBP3-32322] - Western blot analysis of EXOC2 on different lysates with Rabbit anti-EXOC2 antibody (NBP3-32322) at 1/1,000 dilution.


Lane 1: K-562 cell lysate (15 ug/Lane)
 Lane 2: Jurkat cell lysate (15 ug/Lane)
 Lane 3: HepG2 cell lysate (15 ug/Lane)
 Lane 4: MCF7 cell lysate (15 ug/Lane)
 Lane 5: MDA-MB-231 cell lysate (15 ug/Lane)
 Lane 6: HeLa cell lysate (15 ug/Lane)
 Lane 7: HEK-293 cell lysate (15 ug/Lane)
 Lane 8: U-2 OS cell lysate (15 ug/Lane)
 Lane 9: Human brain tissue lysate (30 ug/Lane)
 Lane 10: Mouse brain tissue lysate (30 ug/Lane)
 Lane 11: Rat brain tissue lysate (30 ug/Lane)

Predicted band size: 104 kDa
 Observed band size: 104 kDa

Exposure time: 3 minutes;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-32322) at 1/1,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:100,000 dilution was used for 1 hour at room temperature.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-32322

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-83786PEP	EXOC2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-32322

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

