

# Product Datasheet

## PRKCSH Antibody (A6H9-R)

### NBP3-32386

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-32386](http://www.novusbio.com/NBP3-32386)

Updated 8/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-32386](http://www.novusbio.com/reviews/destination/NBP3-32386)



**NBP3-32386**

PRKCSH Antibody (A6H9-R)

| Product Information         |                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                 |
| Concentration               | 1 mg/ml                                                                                |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality                   | Monoclonal                                                                             |
| Clone                       | A6H9-R                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                     |
| Isotype                     | IgG1                                                                                   |
| Purity                      | Protein A purified                                                                     |
| Buffer                      | PBS (pH7.4), 0.1% BSA and 40% Glycerol                                                 |
| Target Molecular Weight     | 59 kDa                                                                                 |
| Product Description         |                                                                                        |
| Host                        | Mouse                                                                                  |
| Gene ID                     | 5589                                                                                   |
| Gene Symbol                 | PRKCSH                                                                                 |
| Species                     | Human                                                                                  |
| Immunogen                   | Recombinant protein within human PRKCSH aa 51-250. (Uniprot: P14314)                   |
| Product Application Details |                                                                                        |
| Applications                | Western Blot                                                                           |
| Recommended Dilutions       | Western Blot 1:1000                                                                    |



## Images

Western Blot: PRKCSH Antibody (A6H9-R) [NBP3-32386] - Western blot analysis of PRKCSH on different lysates with Mouse anti-PRKCSH antibody (NBP3-32386) at 1/1,000 dilution.

Lane 1: HeLa cell lysate  
 Lane 2: 293T cell lysate  
 Lane 3: Daudi cell lysate  
 Lane 4: Jurkat cell lysate  
 Lane 5: A431 cell lysate


Lysates/proteins at 20 ug/Lane.

Predicted band size: 59 kDa  
 Observed band size: 80 kDa

Exposure time: 5 minutes;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDm/TBST for 1 hour at room temperature. The primary antibody (NBP3-32386) at 1/1,000 dilution was used in 5% NFDm/TBST at 4 overnight. Goat Anti-Mouse IgG - HRP Secondary Antibody at 1/50,000 dilution was used for 1 hour at room temperature.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-32386**

---

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97005-0.5mg   | Mouse IgG1 Isotype Control (MG1)                        |
| H00005589-P02-10ug | Recombinant Human PRKCSH GST (N-Term) Protein           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-32386](http://www.novusbio.com/reviews/submit/NBP3-32386)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

