

Product Datasheet

N-Acetylmannosamine Kinase/GNE Antibody (JE64-81) NBP3-32397

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-32397

Updated 8/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-32397



NBP3-32397**N-Acetylmannosamine Kinase/GNE Antibody (JE64-81)**

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	JE64-81
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	1*TBS (pH7.4), 0.05% BSA and 40% Glycerol
Target Molecular Weight	79 kDa

Product Description	
Host	Rabbit
Gene ID	10020
Gene Symbol	GNE
Species	Human
Immunogen	Synthetic peptide within human N-Acetylmannosamine Kinase/GNE aa 673-722/722. (Uniprot: Q9Y223)

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500



Images

Western Blot: N-Acetylmannosamine Kinase/GNE Antibody (JE64-81) [NBP3-32397] - Western blot analysis of N-Acetylmannosamine Kinase/GNE on different lysates with Rabbit anti-N-Acetylmannosamine Kinase/GNE antibody (NBP3-32397) at 1/500 dilution.

Lane 1: SW480 cell lysate

Lane 2: A549 cell lysate

Lane 3: Hela cell lysate

Lysates/proteins at 10 ug/Lane.


Predicted band size: 79 kDa

Observed band size: 79 kDa

Exposure time: 2 minutes;

10% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-32397) at 1/500 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:300,000 dilution was used for 1 hour at room temperature.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-32397

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-81621PEP	N-Acetylmannosamine Kinase/GNE Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-32397

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

