

# Product Datasheet

## HSP90 beta Antibody (2-1-G3-R)

### NBP3-32446

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-32446](http://www.novusbio.com/NBP3-32446)

Updated 8/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-32446](http://www.novusbio.com/reviews/destination/NBP3-32446)



**NBP3-32446**

HSP90 beta Antibody (2-1-G3-R)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2-1-G3-R
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH7.4), 0.1% BSA and 40% Glycerol
Target Molecular Weight	83.3 kDa

Product Description	
Host	Mouse
Gene ID	3326
Gene Symbol	HSP90AB1
Species	Human, Mouse, Rat
Immunogen	Synthetic peptide within Human HSP90 beta aa 311-360 / 724. (Uniprot: P08238)

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000-1:2000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:5000



## Images

Western Blot: HSP90 beta Antibody (2-1-G3-R) [NBP3-32446] - Western blot analysis of HSP90 beta on different lysates with Mouse anti-HSP90 beta antibody (NBP3-32446) at 1/1,000 dilution.

Lane 1: HeLa cell lysate (20 ug/Lane)  
 Lane 2: HepG2 cell lysate (20 ug/Lane)  
 Lane 3: A549 cell lysate (20 ug/Lane)  
 Lane 4: K-562 cell lysate (20 ug/Lane)  
 Lane 5: Jurkat cell lysate (20 ug/Lane)  
 Lane 6: NIH/3T3 cell lysate (20 ug/Lane)  
 Lane 7: C6 cell lysate (20 ug/Lane)  
 Lane 8: PC-12 cell lysate (20 ug/Lane)  
 Lane 9: Mouse brain tissue lysate (40 ug/Lane)  
 Lane 10: Mouse heart tissue lysate (40 ug/Lane)  
 Lane 11: Rat brain tissue lysate (40 ug/Lane)

Predicted band size: 83 kDa  
 Observed band size: 83 kDa



Exposure time: 53 seconds;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-32446) at 1/1,000 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Mouse IgG - HRP Secondary Antibody at 1/100,000 dilution was used for 1 hour at room temperature.

Immunohistochemistry: HSP90 beta Antibody (2-1-G3-R) [NBP3-32446] - Immunohistochemical analysis of paraffin-embedded human brain tissue with Mouse anti-HSP90 beta antibody (NBP3-32446) at 1/5,000 dilution.


The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH<sub>2</sub>O and PBS, and then probed with the primary antibody (NBP3-32446) at 1/5,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.



Immunocytochemistry/ Immunofluorescence: HSP90 beta Antibody (2-1-G3-R) [NBP3-32446] - Immunocytochemistry analysis of HeLa cells labeling HSP90 beta with Mouse anti-HSP90 beta antibody (NBP3-32446) at 1/100 dilution.

Cells were fixed in 100% precooled methanol for 5 minutes at room temperature, then blocked with 1% BSA in 10% negative goat serum for 1 hour at room temperature. Cells were then incubated with Mouse anti-HSP90 beta antibody (NBP3-32446) at 1/100 dilution in 1% BSA in PBST overnight at 4 °C. Goat Anti-Mouse IgG H&L (iFluor™ 488) was used as the secondary antibody at 1/1,000 dilution. PBS instead of the primary antibody was used as the secondary antibody only control. Nuclear DNA was labelled in blue with DAPI.

beta Tubulin (red) was stained at 1/100 dilution overnight at +4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 594) were used as the secondary antibody at 1/1,000 dilution.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-32446**

---

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-77561PEP	HSP90 beta Antibody Blocking Peptide

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-32446](http://www.novusbio.com/reviews/submit/NBP3-32446)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

