

Product Datasheet

MDH1 Antibody (JE35-03) NBP3-32582

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-32582

Updated 8/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-32582



NBP3-32582

MDH1 Antibody (JE35-03)


Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	JE35-03
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	1*TBS (pH7.4), 0.05% BSA and 40% Glycerol
Target Molecular Weight	36 kDa

Product Description	
Host	Rabbit
Gene ID	4190
Gene Symbol	MDH1
Species	Human, Mouse, Rat
Immunogen	Synthetic peptide. (Uniprot: P40925)


Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Knockout Validated
Recommended Dilutions	Western Blot 1:1000, Flow Cytometry 1:1000, Immunocytochemistry/ Immunofluorescence 1:100, Immunoprecipitation Validated for Immunoprecipitation from YCharOS Inc. (YCharOS.com), Knockout Validated Knockout Validated from YCharOS Inc. (YCharOS.com)

Images


Western blot shows lysates of HAP1 parental cell line and MDH1 knockout HAP1 cell line (KO). Nitrocellulose membrane was probed with MDH1 Antibody (JE35-03) (Catalog # NBP3-32582) followed by HRP-conjugated secondary antibody. A specific band was detected for MDH1 at approximately 36.4 kDa (as indicated) in the parental HAP1 cell line, but is not detectable in knockout HAP1 cell line. Primary antibody dilution used: 1/1000. The Ponceau stained transfer of the blot is shown. This experiment was conducted under reducing conditions. Image, protocol, and testing courtesy of YCharOS Inc. See ycharos.com for additional details.



HAP1 cell line lysates were prepared and immunoprecipitation was performed using 2.0 µg of MDH1 Antibody (JE35-03) (Catalog # NBP3-32582) pre-coupled to Dynabeads Protein A. Immunoprecipitated MDH1 was detected in Western Blot with a rabbit MDH1 antibody used at 1/10000. The Ponceau stained transfer of the blot is shown. SM=4% starting material; UB=4% unbound fraction; IP=immunoprecipitate; HC=antibody heavy chain. Image, protocol and testing courtesy of YCharOS Inc. (ycharos.com).



HAP1 parental cell line WT and MDH1 HAP1 KO cells were labelled with a green or a far-red fluorescent dye, respectively. Cells were stained with MDH1 Antibody (JE35-03) (Catalog # NBP3-32582) followed by incubation with an Alexa-fluor 555 conjugated secondary antibody (upper panel). DAPI-only counterstained cells shown on a lower panel. Acquisition of the blue (nucleus-DAPI), green (identification of WT cells), red (antibody staining) and far-red (identification of KO cells) channels was performed. Representative images of the blue and red (grayscale) channels are shown. WT and KO cells are outlined with green and magenta dashed line, respectively. Primary antibody dilution used: 1/500. Image, protocol and testing courtesy of YCharOS Inc. (ycharos.com).



Western Blot: MDH1 Antibody (JE35-03) [NBP3-32582] - Western blot analysis of MDH1 on different lysates with Rabbit anti-MDH1 antibody (NBP3-32582) at 1/1,000 dilution.


Lane 1: Jurkat cell lysate (20 ug/Lane)
Lane 2: HepG2 cell lysate (20 ug/Lane)
Lane 3: U-937 cell lysate (20 ug/Lane)
Lane 4: HL-60 cell lysate (20 ug/Lane)
Lane 5: NIH/3T3 cell lysate (20 ug/Lane)
Lane 6: Rat kidney tissue lysate (40 ug/Lane)

Predicted band size: 36 kDa
Observed band size: 35 kDa

Exposure time: 6 seconds; ECL;


4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-32582) at 1/1,000 dilution was used in 5% NFDM/TBST at 4 overnight. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1/50,000 dilution was used for 1 hour at room temperature.



Immunocytochemistry/ Immunofluorescence: MDH1 Antibody (JE35-03) [NBP3-32582] - Immunocytochemistry analysis of U-937 cells labeling MDH1 with Rabbit anti-MDH1 antibody (NBP3-32582) at 1/100 dilution.


Cells were fixed in 4% paraformaldehyde for 20 minutes at room temperature, permeabilized with 0.1% Triton X-100 in PBS for 5 minutes at room temperature, then blocked with 1% BSA in 10% negative goat serum for 1 hour at room temperature. Cells were then incubated with Rabbit anti-MDH1 antibody (NBP3-32582) at 1/100 dilution in 1% BSA in PBST overnight at 4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 488) was used as the secondary antibody at 1/1,000 dilution. PBS instead of the primary antibody was used as the secondary antibody only control. Nuclear DNA was labelled in blue with DAPI.



Beta tubulin (red) was stained at 1/100 dilution overnight at +4 °C. Goat Anti-Mouse IgG H&L (iFluor™ 594) was used as the secondary antibody at 1/1,000 dilution.

Flow Cytometry: MDH1 Antibody (JE35-03) [NBP3-32582] - Flow cytometric analysis of U-937 cells labeling MDH1.

Cells were fixed and permeabilized. Then stained with the primary antibody (NBP3-32582, 1 µg/mL) (red) compared with Rabbit IgG Isotype Control (green). After incubation of the primary antibody at +4 °C for an hour, the cells were stained with a iFluor™ 488 conjugate-Goat anti-Rabbit IgG Secondary antibody at 1/1,000 dilution for 30 minutes at +4 °C. Unlabelled sample was used as a control (cells without incubation with primary antibody; black).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-32582

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-45315-0.1mg	Recombinant Human MDH1 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-32582

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

