

# Product Datasheet

## SSX2IP Antibody (JE65-50)

### NBP3-32978

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-32978](http://www.novusbio.com/NBP3-32978)

Updated 8/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-32978](http://www.novusbio.com/reviews/destination/NBP3-32978)



**NBP3-32978**

SSX2IP Antibody (JE65-50)

| Product Information     |                                                                                        |
|-------------------------|----------------------------------------------------------------------------------------|
| Unit Size               | 100 ul                                                                                 |
| Concentration           | 1 mg/ml                                                                                |
| Storage                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality               | Monoclonal                                                                             |
| Clone                   | JE65-50                                                                                |
| Preservative            | 0.05% Sodium Azide                                                                     |
| Isotype                 | IgG                                                                                    |
| Purity                  | Protein A purified                                                                     |
| Buffer                  | 1*TBS (pH7.4), 0.05% BSA and 40% Glycerol                                              |
| Target Molecular Weight | 71 kDa                                                                                 |

| Product Description |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Host                | Rabbit                                                                    |
| Gene ID             | 117178                                                                    |
| Gene Symbol         | SSX2IP                                                                    |
| Species             | Human                                                                     |
| Immunogen           | Synthetic peptide within human SLC25A12 aa 565-614/614. (Uniprot: Q9Y2D8) |

| Product Application Details |                    |
|-----------------------------|--------------------|
| Applications                | Western Blot       |
| Recommended Dilutions       | Western Blot 1:500 |

**Images**

Western Blot: SSX2IP Antibody (JE65-50) [NBP3-32978] - Western blot analysis of SSX2IP on different lysates with Rabbit anti-SSX2IP antibody (NBP3-32978) at 1/500 dilution.

Lane 1: Daudi cell lysate  
Lane 2: K562 cell lysate


Lysates/proteins at 10 ug/Lane.

Predicted band size: 71 kDa  
Observed band size: 80/90 kDa

Exposure time: 2 minutes;

10% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-32978) at 1/500 dilution was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:300,000 dilution was used for 1 hour at room temperature.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-32978**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-85956PEP | SSX2IP Recombinant Protein Antigen                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-32978](http://www.novusbio.com/reviews/submit/NBP3-32978)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

