

Product Datasheet

UQCRC1 Antibody (PSH01-23)

NBP3-33078

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-33078

Updated 8/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-33078



NBP3-33078

UQCRC1 Antibody (PSH01-23)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PSH01-23
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH7.4), 0.1% BSA and 40% Glycerol
Target Molecular Weight	53 kDa
Product Description	
Host	Rabbit
Gene ID	7384
Gene Symbol	UQCRC1
Species	Human, Mouse, Rat
Immunogen	Recombinant protein within human UQCRC1 aa 1-480 / 480. (Uniprot: P31930)
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:5000, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 1:100, Immunohistochemistry-Paraffin 1:5000



Images

Western Blot: UQCRC1 Antibody (PSH01-23) [NBP3-33078] - Western blot analysis of UQCRC1 on different lysates with Rabbit anti-UQCRC1 antibody (NBP3-33078) at 1/5,000 dilution.

Lane 1: HeLa cell lysate (20 ug/Lane)
 Lane 2: HepG2 cell lysate (20 ug/Lane)
 Lane 3: HCT 116 cell lysate (20 ug/Lane)
 Lane 4: 293T cell lysate (20 ug/Lane)
 Lane 5: PC-3M cell lysate (20 ug/Lane)
 Lane 6: Jurkat cell lysate (20 ug/Lane)
 Lane 7: NIH/3T3 cell lysate (20 ug/Lane)
 Lane 8: PC-12 cell lysate (20 ug/Lane)
 Lane 9: Mouse kidney tissue lysate (40 ug/Lane)
 Lane 10: Mouse brain tissue lysate (40 ug/Lane)
 Lane 11: Rat kidney tissue lysate (40 ug/Lane)
 Lane 12: Rat brain tissue lysate (40 ug/Lane)

Predicted band size: 53 kDa
 Observed band size: 45 kDa

Exposure time: 8 seconds;

4-20% SDS-PAGE gel.




Proteins were transferred to a PVDF membrane and blocked with 5% NFDm/TBST for 1 hour at room temperature. The primary antibody (NBP3-33078) at 1/5,000 dilution was used in 5% NFDm/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1:100,000 dilution was used for 1 hour at room temperature.

Immunohistochemistry: UQCRC1 Antibody (PSH01-23) [NBP3-33078] - Immunohistochemical analysis of paraffin-embedded human kidney tissue with Rabbit anti-UQCRC1 antibody (NBP3-33078) at 1/5,000 dilution.

The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 1% BSA for 20 minutes at room temperature, washed with ddH₂O and PBS, and then probed with the primary antibody (NBP3-33078) at 1/5,000 dilution for 1 hour at room temperature. The detection was performed using an HRP conjugated compact polymer system. DAB was used as the chromogen. Tissues were counterstained with hematoxylin and mounted with DPX.

Immunohistochemistry: UQCRC1 Antibody (PSH01-23) [NBP3-33078] - Immunofluorescence analysis of paraffin-embedded human liver tissue labeling UQCRC1 with Rabbit anti-UQCRC1 antibody (NBP3-33078) at 1/1,000 dilution.


The section was pre-treated using heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 20 minutes. The tissues were blocked in 10% negative goat serum for 1 hour at room temperature, washed with PBS, and then probed with the primary antibody (NBP3-33078, green) at 1/1,000 dilution overnight at 4 °C, washed with PBS. Goat Anti-Rabbit IgG H&L (iFluor™ 488) was used as the secondary antibody at 1/1,000 dilution. Nuclei were counterstained with DAPI (blue).



Immunocytochemistry/ Immunofluorescence: UQCRC1 Antibody (PSH01-23) [NBP3-33078] - Immunocytochemistry analysis of HepG2 cells labeling UQCRC1 with Rabbit anti-UQCRC1 antibody (NBP3-33078) at 1/100 dilution.

Cells were fixed in 4% paraformaldehyde for 20 minutes at room temperature, permeabilized with 0.1% Triton X-100 in PBS for 5 minutes at room temperature, then blocked with 1% BSA in 10% negative goat serum for 1 hour at room temperature. Cells were then incubated with Rabbit anti-UQCRC1 antibody (NBP3-33078) at 1/100 dilution in 1% BSA in PBST overnight at 4 °C. Goat Anti-Rabbit IgG H&L (iFluor™ 488) was used as the secondary antibody at 1/1,000 dilution. PBS instead of the primary antibody was used as the secondary antibody only control. Nuclear DNA was labelled in blue with DAPI.

Beta tubulin (red) was stained at 1/100 dilution overnight at +4 °C. Goat Anti-Mouse IgG H&L (iFluor™ 594) was used as the secondary antibody at 1/1,000 dilution.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-33078

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-51608-0.1mg	Recombinant Human UQCRC1 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-33078

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

