

Product Datasheet

p38 beta/gamma/delta Antibody (PSH03-59) **NBP3-33145**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-33145

Updated 8/7/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-33145



NBP3-33145

p38 beta/gamma/delta Antibody (PSH03-59)

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PSH03-59
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH7.4), 0.1% BSA and 40% Glycerol
Target Molecular Weight	42 kDa

Product Description	
Host	Rabbit
Gene ID	5600
Gene Symbol	MAPK11
Species	Human, Mouse, Rat, Monkey
Specificity/Sensitivity	Anti-p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Antibody [PSH03-59] detects endogenous levels of total p38 beta, p38 gamma and p38 delta MAPK protein. This antibody does not recognize p38 alpha MAPK protein.
Immunogen	Recombinant protein within human MAPK12 aa 1-367 / 367. (Uniprot: Q15759)

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:2000



Images

Western Blot: p38 beta/gamma/delta Antibody (PSH03-59) [NBP3-33145] - Western blot analysis of p38 beta/gamma/delta on different lysates with Rabbit anti-p38 beta/gamma/delta antibody (NBP3-33145) at 1/2,000 dilution.

Lane 1: 293T cell lysate
 Lane 2: HeLa cell lysate
 Lane 3: MCF7 cell lysate
 Lane 4: K-562 cell lysate
 Lane 5: MDA-MB-231 cell lysate
 Lane 6: COS-1 cell lysate
 Lane 7: RAW264.7 cell lysate
 Lane 8: Neuro-2a cell lysate
 Lane 9: PC-12 cell lysate


Lysates/proteins at 30 ug/Lane.

Predicted band size: 42 kDa
 Observed band size: 42 kDa

Exposure time: 1 minute 2 seconds;

4-20% SDS-PAGE gel.

Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (NBP3-33145) at 1/2,000 dilution was used in 5% NFDM/TBST at 4C overnight. Goat Anti-Rabbit IgG - HRP Secondary Antibody at 1/50,000 dilution was used for 1 hour at room temperature.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-33145

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-33145

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

