

# Product Datasheet

## IKBKAP Antibody [PE] NBP3-35099PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-35099PE](http://www.novusbio.com/NBP3-35099PE)

Updated 3/24/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-35099PE](http://www.novusbio.com/reviews/destination/NBP3-35099PE)



**NBP3-35099PE**

IKBKAP Antibody [PE]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                   |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                         |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Conjugate                   | PE                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Purity                      | Affinity purified                                                                                                                                                                                                                                                                                                                                                                                                          |
| Buffer                      | PBS                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Gene ID                     | 8518                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Gene Symbol                 | ELP1                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                          |
| Immunogen                   | <p>Recombinant fusion protein containing a sequence corresponding to amino acids 1093-1332 of human IKBKAP (NP_003631.2).</p> <p>Sequence:<br/> LRLVYKYNRLDIIETNVKPSILEAQKNYMAFLDSQTATFSRHKKRLLVVRELKEQ<br/> AQQAGLDDEVPHGQESDLFSETSSVSGSEMSGKYSHSNSRISARSSKNRRK<br/> AERKKHSLKEGSPLEDLALLEALSEVVQNTENLKDEVYHILKVLFLFEFDEQGR<br/> ELQKAFEDTLQLMERSLPEIWTLYQQNSATPVLGPNSTANSIMASYQQQKTS<br/> VPVLDAELFIPPKINRRTQWKLSLLD</p> |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Applications                | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                   |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                   |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-35099PE**

---

|               |                                    |
|---------------|------------------------------------|
| NBP2-24983    | Rabbit IgG Isotype Control [PE]    |
| NBP2-56082PEP | IKBKAP Recombinant Protein Antigen |
| 210-TA-005    | TNF-alpha [Unconjugated]           |
| NBP2-07968    | IKBKAP Overexpression Lysate       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-35099PE](http://www.novusbio.com/reviews/submit/NBP3-35099PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

