

# Product Datasheet

## Wnt-11 Antibody (D13) NBP3-42179

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-42179](http://www.novusbio.com/NBP3-42179)

Updated 11/11/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-42179](http://www.novusbio.com/reviews/destination/NBP3-42179)



**NBP3-42179**

Wnt-11 Antibody (D13)

| Product Information         |                                                                                                            |
|-----------------------------|------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                     |
| Concentration               | 1 mg/ml                                                                                                    |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                     |
| Clonality                   | Monoclonal                                                                                                 |
| Clone                       | D13                                                                                                        |
| Preservative                | 0.05% Proclin 300                                                                                          |
| Isotype                     | IgG2b Kappa                                                                                                |
| Purity                      | Protein A or G purified                                                                                    |
| Buffer                      | 0.01M PBS, pH7.4, 50% glycerol                                                                             |
| Product Description         |                                                                                                            |
| Host                        | Mouse                                                                                                      |
| Gene ID                     | 7481                                                                                                       |
| Gene Symbol                 | WNT11                                                                                                      |
| Species                     | Human, Porcine                                                                                             |
| Immunogen                   | Immunogen sequence: Ser204-Glu349 (Uniprot#: O96014)                                                       |
| Product Application Details |                                                                                                            |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                |
| Recommended Dilutions       | Western Blot 0.5-2ug/mL, Immunohistochemistry 5-20ug/mL, Immunocytochemistry/ Immunofluorescence 5-20ug/mL |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-42179**

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-43317-0.5mg | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)    |
| NBP2-13523PEP    | Wnt-11 Recombinant Protein Antigen                      |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-42179](http://www.novusbio.com/reviews/submit/NBP3-42179)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

