

# Product Datasheet

## Cylicin 1 Recombinant Protein Antigen NBP3-43686PEP

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-43686PEP](http://www.novusbio.com/NBP3-43686PEP)

Updated 2/3/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-43686PEP](http://www.novusbio.com/reviews/destination/NBP3-43686PEP)



**NBP3-43686PEP****Cylicin 1 Recombinant Protein Antigen**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 100 ul                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Preservative</b>                | No Preservative                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>                      | >80% by SDS-PAGE and Coomassie blue staining                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Buffer</b>                      | PBS and 1M Urea, pH 7.4.                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Description</b>                 | <p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human Cylicin 1</p> <p><b>Source:</b> <i>E.coli</i></p> <p><b>Amino Acid Sequence:</b><br/>GAKKKIDESDGTSA NSKMEGLESKRGFRMSSKKTTTFNEKG EKASTGRVPPSR<br/>EKPPLPACEPSLPSPKVRRLCWCKMPPPPPKPRYAP</p> <p><b>Fusion Tag:</b> N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p>                                                           |
| <b>Gene ID</b>                     | 1538                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Gene Symbol</b>                 | CYLC1                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Species</b>                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Applications</b>                | Antibody Competition                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Recommended Dilutions</b>       | Antibody Competition 10-100 molar excess                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Application Notes</b>           | <p>This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP3-43686It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml.For current lot information, including availability, please contact our technical support team click <a href="mailto:nb-technical@bio-techne.com">nb-technical@bio-techne.com</a>.</p> |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-43686PEP](http://www.novusbio.com/reviews/submit/NBP3-43686PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

