

# Product Datasheet

## MLH3 Antibody (CL14084) [FITC]

### NBP3-43849F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-43849F](http://www.novusbio.com/NBP3-43849F)

Updated 3/27/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-43849F](http://www.novusbio.com/reviews/destination/NBP3-43849F)



**NBP3-43849F**

MLH3 Antibody (CL14084) [FITC]

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                                         |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | CL14084                                                                                                          |
| Preservative        | 0.05% Sodium Azide                                                                                               |
| Isotype             | IgG2b                                                                                                            |
| Conjugate           | FITC                                                                                                             |
| Purity              | Protein A purified                                                                                               |
| Buffer              | PBS                                                                                                              |

| Product Description |                                                                                                                           |
|---------------------|---------------------------------------------------------------------------------------------------------------------------|
| Host                | Mouse                                                                                                                     |
| Gene ID             | 27030                                                                                                                     |
| Gene Symbol         | MLH3                                                                                                                      |
| Species             | Human                                                                                                                     |
| Immunogen           | This antibody was generated using a synthetic peptide of Q9UHC1, with the exact immunogen sequence remaining proprietary. |

| Product Application Details |                                                                                      |
|-----------------------------|--------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.               |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-43849F**

---

|               |                                                    |
|---------------|----------------------------------------------------|
| NBP2-27229    | Mouse IgG2b Isotype Control (MPC-11) [FITC]        |
| NBP2-55958PEP | MLH3 Recombinant Protein Antigen                   |
| NBP2-07349    | MLH3 Overexpression Lysate                         |
| AF2535        | MMR/CD206/Mannose Receptor Antibody [Unconjugated] |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-43849F](http://www.novusbio.com/reviews/submit/NBP3-43849F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

