

# Product Datasheet

## LONRF3 Antibody (OTI7F7)

### NBP3-44961-100ul

Unit Size: 100 ul

Store at -20 degrees C. Avoid freeze/thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-44961](http://www.novusbio.com/NBP3-44961)

Updated 3/6/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-44961](http://www.novusbio.com/reviews/destination/NBP3-44961)



**NBP3-44961-100ul**

LONRF3 Antibody (OTI7F7)

| Product Information         |                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                |
| Storage                     | Store at -20 degrees C. Avoid freeze/thaw cycles.                                                                          |
| Clonality                   | Monoclonal                                                                                                                 |
| Clone                       | OTI7F7                                                                                                                     |
| Preservative                | 0.02% Sodium Azide                                                                                                         |
| Isotype                     | IgG1                                                                                                                       |
| Purity                      | >95%                                                                                                                       |
| Buffer                      | PBS (pH 7.3) 1% BSA, 50% glycerol                                                                                          |
| Product Description         |                                                                                                                            |
| Host                        | Mouse                                                                                                                      |
| Gene ID                     | 79836                                                                                                                      |
| Gene Symbol                 | LONRF3                                                                                                                     |
| Species                     | Human                                                                                                                      |
| Immunogen                   | Human recombinant protein fragment corresponding to amino acids 206-456 of human LONRF3 (NP_001027026) produced in E.coli. |
| Product Application Details |                                                                                                                            |
| Applications                | Western Blot                                                                                                               |
| Recommended Dilutions       | Western Blot 1:500                                                                                                         |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-44961-100ul**

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97005-0.5mg | Mouse IgG1 Isotype Control (MG1)                        |
| NBP2-68769PEP    | LONRF3 Recombinant Protein Antigen                      |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-44961](http://www.novusbio.com/reviews/submit/NBP3-44961)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

