

# Product Datasheet

## L1TD1 Antibody (L1TD1/7941)

### NBP3-48954-100ug

Unit Size: 100 ug

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-48954](http://www.novusbio.com/NBP3-48954)

Updated 8/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-48954](http://www.novusbio.com/reviews/destination/NBP3-48954)



**NBP3-48954-100ug**

L1TD1 Antibody (L1TD1/7941)

| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 100 ug                                                                                      |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C.                                                                                |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | L1TD1/7941                                                                                  |
| Preservative            | 0.05% Sodium Azide                                                                          |
| Isotype                 | IgG2 Kappa                                                                                  |
| Purity                  | Protein A purified                                                                          |
| Buffer                  | 10 mM PBS with 0.05% BSA                                                                    |
| Target Molecular Weight | 99 kDa                                                                                      |

| Product Description |                                                                                                                                                                                                                                                                                                                                           |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | 200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.<br><br>Antibody with azide - store at 2 to 8 C. Antibody without azide - store at -20 to -80 C. Antibody is stable for 24 months. Non-hazardous. No MSDS required. |
| Host                | Mouse                                                                                                                                                                                                                                                                                                                                     |
| Gene ID             | 54596                                                                                                                                                                                                                                                                                                                                     |
| Gene Symbol         | L1TD1                                                                                                                                                                                                                                                                                                                                     |
| Species             | Human                                                                                                                                                                                                                                                                                                                                     |
| Immunogen           | Recombinant fragment of human L1TD1 protein (exact sequence is proprietary) (Uniprot: Q5T7N2)                                                                                                                                                                                                                                             |

| Product Application Details |                      |
|-----------------------------|----------------------|
| Applications                | Immunohistochemistry |
| Recommended Dilutions       | Immunohistochemistry |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-48954-100ug**

---

|            |                                                         |
|------------|---------------------------------------------------------|
| HAF007     | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B    | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| 314-BP-010 | BMP-4 [Unconjugated]                                    |
| AF2018     | SOX2 Antibody [Unconjugated]                            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-48954](http://www.novusbio.com/reviews/submit/NBP3-48954)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

