

# Product Datasheet

## EGFR [p Tyr1173] Antibody (24GB4400) NBP3-50794-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-50794](http://www.novusbio.com/NBP3-50794)

Updated 3/24/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-50794](http://www.novusbio.com/reviews/destination/NBP3-50794)



**NBP3-50794-100ul**

EGFR [p Tyr1173] Antibody (24GB4400)

| Product Information |   |
|---------------------|---|
| Unit Size           | 100 ul  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at -20C. Avoid freeze-thaw cycles.  |
| Clonality           | Monoclonal  |
| Clone               | 24GB4400  |
| Preservative        | 0.02% Sodium Azide  |
| Isotype             | IgG   |
| Purity              | Affinity purified   |
| Buffer              | PBS (pH 7.4), 50% glycerol  |

| Product Description     |   |
|-------------------------|---|
| Host                    | Rabbit  |
| Gene ID                 | 1956  |
| Gene Symbol             | EGFR  |
| Species                 | Human, Rat  |
| Specificity/Sensitivity | This antibody detects endogenous levels of EGFR only when phosphorylated at Y1173 |
| Immunogen               | A synthesized peptide derived from human Phospho-EGFR (Tyr1173) (Uniprot: P00533) |

| Product Application Details |                              |
|-----------------------------|------------------------------|
| Applications                | Western Blot, Flow Cytometry |
| Recommended Dilutions       | Western Blot, Flow Cytometry |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-50794-100ul**

---

|               |   |
|---------------|---|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-84815PEP | EGFR Recombinant Protein Antigen                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-50794](http://www.novusbio.com/reviews/submit/NBP3-50794)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

