

# Product Datasheet

## ErbB2/Her2 Antibody (24GB5655) NBP3-51733-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-51733](http://www.novusbio.com/NBP3-51733)

Updated 3/24/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-51733](http://www.novusbio.com/reviews/destination/NBP3-51733)



**NBP3-51733-100ul**

ErbB2/Her2 Antibody (24GB5655)

| Product Information         |                                                                                                                                                                                                                              |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                                                       |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                  |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                     |
| Clonality                   | Monoclonal                                                                                                                                                                                                                   |
| Clone                       | 24GB5655                                                                                                                                                                                                                     |
| Preservative                | 0.02% Sodium Azide                                                                                                                                                                                                           |
| Isotype                     | IgG2a                                                                                                                                                                                                                        |
| Purity                      | Affinity purified                                                                                                                                                                                                            |
| Buffer                      | PBS (pH 7.4), 50% glycerol                                                                                                                                                                                                   |
| Product Description         |                                                                                                                                                                                                                              |
| Description                 | Reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Host                        | Mouse                                                                                                                                                                                                                        |
| Gene ID                     | 2064                                                                                                                                                                                                                         |
| Gene Symbol                 | ERBB2                                                                                                                                                                                                                        |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                            |
| Immunogen                   | Recombinant protein of human ErbB2/Her2 (Uniprot: P04626)                                                                                                                                                                    |
| Product Application Details |                                                                                                                                                                                                                              |
| Applications                | Western Blot, Flow Cytometry, Knockdown Validated                                                                                                                                                                            |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Knockdown Validated                                                                                                                                                                            |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-51733-100ul**

|               |                                                    |
|---------------|----------------------------------------------------|
| HAF007        | Goat anti-Mouse IgG Secondary Antibody [HRP]       |
| NB7539        | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-96778    | Mouse IgG2a Isotype Control (M2A)                  |
| NBP2-29624PEP | ErbB2/Her2 Antibody Blocking Peptide               |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-51733](http://www.novusbio.com/reviews/submit/NBP3-51733)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

