

# Product Datasheet

## MIG2/Kindlin-2 Antibody (24GB6190) NBP3-52047-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-52047](http://www.novusbio.com/NBP3-52047)

Updated 3/24/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-52047](http://www.novusbio.com/reviews/destination/NBP3-52047)



**NBP3-52047-100ul**

MIG2/Kindlin-2 Antibody (24GB6190)

| Product Information         |                                                                                                                                                                                                                              |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                                                       |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                  |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                     |
| Clonality                   | Monoclonal                                                                                                                                                                                                                   |
| Clone                       | 24GB6190                                                                                                                                                                                                                     |
| Preservative                | 0.02% Sodium Azide                                                                                                                                                                                                           |
| Isotype                     | IgG2a                                                                                                                                                                                                                        |
| Purity                      | Affinity purified                                                                                                                                                                                                            |
| Buffer                      | PBS (pH 7.4), 50% glycerol                                                                                                                                                                                                   |
| Product Description         |                                                                                                                                                                                                                              |
| Description                 | Reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Host                        | Mouse                                                                                                                                                                                                                        |
| Gene ID                     | 10979                                                                                                                                                                                                                        |
| Gene Symbol                 | FERMT2                                                                                                                                                                                                                       |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                            |
| Immunogen                   | Recombinant protein of human MIG2/Kindlin-2 (Uniprot: Q96AC1)                                                                                                                                                                |
| Product Application Details |                                                                                                                                                                                                                              |
| Applications                | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Knockout Validated                                                                                                                                    |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Knockout Validated                                                                                                                                    |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-52047-100ul**

|                    |                                                       |
|--------------------|-------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]          |
| NB7539             | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]    |
| NBP1-96778         | Mouse IgG2a Isotype Control (M2A)                     |
| H00010979-P01-10ug | Recombinant Human MIG2/Kindlin-2 GST (N-Term) Protein |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-52047](http://www.novusbio.com/reviews/submit/NBP3-52047)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

