

# Product Datasheet

## **LYPD6B Antibody (634CT7.10.1) [mFluor Violet 610 SE] NBP3-52310MFV610**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-52310MFV610](http://www.novusbio.com/NBP3-52310MFV610)

Updated 7/20/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-52310MFV610](http://www.novusbio.com/reviews/destination/NBP3-52310MFV610)



**NBP3-52310MFV610**

LYPD6B Antibody (634CT7.10.1) [mFluor Violet 610 SE]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | This product is unpurified. The exact concentration of antibody is not quantifiable.                                                                                                                                                                                                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | 634CT7.10.1                                                                                                                                                                                                                                                                                                                                     |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                                                                                                           |
| Conjugate                   | mFluor Violet 610 SE                                                                                                                                                                                                                                                                                                                            |
| Purity                      | Ascites                                                                                                                                                                                                                                                                                                                                         |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 130576                                                                                                                                                                                                                                                                                                                                          |
| Gene Symbol                 | LYPD6B                                                                                                                                                                                                                                                                                                                                          |
| Species                     | Human, Mouse                                                                                                                                                                                                                                                                                                                                    |
| Reactivity Notes            | Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions. |
| Immunogen                   | This LYPD6B antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 79-108 amino acids from the Central region of human LYPD6B.(Accession#:NP_808879.2)                                                                                                                                                       |
| Notes                       | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.                                                                                                           |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                    |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                                                    |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-52310MFV610**

|                  |                                                             |
|------------------|-------------------------------------------------------------|
| NBP2-27231MFV610 | Mouse IgG2b Isotype Control (MPC-11) [mFluor Violet 610 SE] |
| DCF300           | Coagulation Factor III/Tissue Factor [HRP]                  |
| AF904            | Contactin-1 Antibody [Unconjugated]                         |
| NBP2-67471       | c-jun [p Ser63] Antibody (SY0297)                           |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-52310MFV610](http://www.novusbio.com/reviews/submit/NBP3-52310MFV610)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

