

# Product Datasheet

## **LLGL1/LLGL2 [p Ser655, p Ser645, p Ser649, p Ser659] Antibody (RB7629) [DyLight 755] NBP3-52327IR**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-52327IR](http://www.novusbio.com/NBP3-52327IR)

Updated 7/20/2026 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-52327IR](http://www.novusbio.com/reviews/destination/NBP3-52327IR)



**NBP3-52327IR**

LLGL1/LLGL2 [p Ser655, p Ser645, p Ser649, p Ser659] Antibody (RB7629) [DyLight 755]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Polyclonal                                                                                  |
| Clone               | RB7629                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Conjugate           | DyLight 755                                                                                 |
| Purity              | Antigen and protein A Affinity-purified                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description     |                                                                                                                                                                                                                                                 |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                    | Rabbit                                                                                                                                                                                                                                          |
| Gene ID                 | 3993                                                                                                                                                                                                                                            |
| Gene Symbol             | LLGL2                                                                                                                                                                                                                                           |
| Species                 | Human                                                                                                                                                                                                                                           |
| Specificity/Sensitivity | this detects phosphorylated Ser 655 and Ser 659 for LLGL1, and Ser 645 and Ser 649 for LLGL2                                                                                                                                                    |
| Immunogen               | This LLGL1/2 Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding S645/S649 of human LLGL1/2.<br>(Accession#:NP_001015002.1;NP_001026973.1;NP_004515.2) |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                               |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot                                                           |
| Recommended Dilutions       | Western Blot                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-52327IR**

---

|              |                                              |
|--------------|----------------------------------------------|
| NBP2-24891IR | Rabbit IgG Isotype Control [DyLight 755]     |
| 236-EG-200   | EGF [Unconjugated]                           |
| AF748        | E-Cadherin Antibody [Unconjugated]           |
| NBP1-85047   | Tight Junction Protein 1 Antibody - BSA Free |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-52327IR](http://www.novusbio.com/reviews/submit/NBP3-52327IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

