

Product Datasheet

CACNA2D2 Antibody (N492/39) - Chimeric - Azide and BSA Free **NBP3-52363**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-52363

Updated 6/11/2026 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-52363



NBP3-52363

CACNA2D2 Antibody (N492/39) - Chimeric - Azide and BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	N492/39
Preservative	No Preservative
Isotype	IgY
Purity	Affinity purified
Buffer	PBS
Target Molecular Weight	130 kDa


Product Description	
Description	This recombinant antibody is a chimeric antibody created by replacing the mouse heavy and light constant regions of clone N492/39 with chicken IgY heavy and light constant regions. As such this antibody retains the same binding performance as the original clone N492/39 but can be detected using standard anti-chicken secondary antibodies allowing flexibility for multiplexing applications. This antibody is expressed recombinantly in HEK293 cells and then affinity purified from the cell culture media and buffer exchanged into 1X PBS.
Host	Chicken
Gene ID	9254
Gene Symbol	CACNA2D2
Species	Human, Mouse, Rat
Immunogen	Recombinant protein amino acids 61-197 (extracellular domain) of mouse CACNA2D2 (Uniprot#: Q6PHS9) produced recombinantly in E.Coli.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin




Images

Representative image of formalin fixed paraffin-embedded section of rat cerebellum section stained with NBP3-52362 (the BSA Free version of this antibody) (brown) and counterstained with hematoxylin (blue). This antibody specifically stains the Purkinje cells in the cerebellum.



Western blotting of rat brain membrane with NBP3-52362 (the BSA Free version of this antibody), recognizing endogenous CCACNA2D2 at ~130 kDa, and detected with anti-chicken-HRP.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-52363

BAF010	Goat anti-Chicken IgY Secondary Antibody [Biotin]
NB7276	Goat anti-Chicken IgM Heavy Chain Secondary Antibody
NBP1-81501PEP	CACNA2D2 Recombinant Protein Antigen
NB100-565	FUS Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-52363

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

