

Product Datasheet

ODZ3 Antibody (N497/10) - Chimeric - Azide and BSA Free NBP3-52367

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-52367

Updated 6/11/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-52367



NBP3-52367


ODZ3 Antibody (N497/10) - Chimeric - Azide and BSA Free

Product Information	
Unit Size	100 ul
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	N497/10
Preservative	No Preservative
Isotype	IgY
Purity	Affinity purified
Buffer	PBS
Target Molecular Weight	303 kDa
Product Description	
Description	This recombinant antibody is a chimeric antibody created by replacing the mouse heavy and light constant regions of clone N497/10 with chicken IgY heavy and light constant regions. As such this antibody retains the same binding performance as the original clone N497/10 but can be detected using standard anti-chicken secondary antibodies allowing flexibility for multiplexing applications. This antibody is expressed recombinantly in HEK293 cells and then affinity purified from the cell culture media and buffer exchanged into 1X PBS.
Host	Chicken
Gene ID	55714
Gene Symbol	TENM3
Species	Human, Mouse, Rat
Immunogen	Recombinant protein comprising amino acids aa 161-201 and 343-396 (extracellular and cytoplasmic domains) of mouse ODZ3 (Uniprot#: Q9WTS6) connected by a Gly8 linker, produced recombinantly in E.Coli.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin




Images

Representative image of formalin fixed paraffin embedded section of rat cerebellum section stained with NBP3-52366 (the BSA Free version of this antibody) (brown) and counterstained with hematoxylin (blue). Magnification: 200X



Western blotting of rat brain membrane with ODZ3 mouse monoclonal antibody (A) and with NBP3-52366 (the BSA Free version of this antibody) (B) recognizing endogenous ODZ3 at ~300 kDa with similar sensitivity. The mouse monoclonal antibody and chicken monoclonal antibody were detected with anti-mouse-HRP and anti-chicken-HRP, respectively.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-52367

BAF010	Goat anti-Chicken IgY Secondary Antibody [Biotin]
NB7276	Goat anti-Chicken IgM Heavy Chain Secondary Antibody
NBP2-13687PEP	ODZ3 Recombinant Protein Antigen
236-EG-200	EGF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-52367

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

