

Product Datasheet

Ankyrin 3 Antibody (N106/36) [FL650] - Chimeric - BSA Free NBP3-52382

Unit Size: 200 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-52382

Updated 6/14/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-52382



NBP3-52382

Ankyrin 3 Antibody (N106/36) [FL650] - Chimeric - BSA Free

Product Information	
Unit Size	200 ul
Concentration	0.5 mg/ml
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	N106/36
Preservative	0.05% Sodium Azide
Isotype	IgY
Conjugate	FL650
Purity	Affinity purified
Buffer	1X PBS (pH 7.4)

Product Description	
Description	<p>FL650 Excitation: 654 nm, Emission: 672nm. Each order is a unique lot. For specialized requirements, please inquire for available options.</p> <p>This recombinant antibody is a chimeric antibody created by replacing the mouse heavy and light constant regions of clone N106/36 with chicken IgY heavy and light constant regions. As such this antibody retains the same binding performance as the original clone N106/36 but can be detected using standard anti-chicken secondary antibodies allowing flexibility for multiplexing applications. This antibody is expressed recombinantly in HEK293 cells and then affinity purified from the cell culture media and conjugated.</p>
Host	Chicken
Gene ID	288
Gene Symbol	ANK3
Species	Human, Mouse, Rat, Feline, Hamster, Primate
Specificity/Sensitivity	Does not cross-react with Ankyrin-B
Immunogen	Fusion protein ~1000 amino acids of Ankyrin 3 (Uniprot#: E9PE32) produced recombinantly in E. Coli

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-52382

H00006326-Q01-10ug	Recombinant Human SCN2A GST (N-Term) Protein
AF3235	Neurofascin Antibody
AF1916	p63/TP73L Antibody [Unconjugated]
NB100-56603	Androgen R/NR3C4 [p Ser213, p Ser210] Antibody (156C135.2) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-52382

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

