

Product Datasheet

ODZ3 Antibody (N497/10) [FL650] - Chimeric - BSA Free NBP3-52434

Unit Size: 200 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-52434

Updated 6/14/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-52434



NBP3-52434

ODZ3 Antibody (N497/10) [FL650] - Chimeric - BSA Free

Product Information	
Unit Size	200 ul
Concentration	0.5 mg/ml
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	N497/10
Preservative	0.05% Sodium Azide
Isotype	IgY
Conjugate	FL650
Purity	Affinity purified
Buffer	1X PBS (pH 7.4)
Target Molecular Weight	303 kDa

Product Description	
Description	<p>FL650 Excitation: 654 nm, Emission: 672nm. Each order is a unique lot. For specialized requirements, please inquire for available options.</p> <p>This recombinant antibody is a chimeric antibody created by replacing the mouse heavy and light constant regions of clone N497/10 with chicken IgY heavy and light constant regions. As such this antibody retains the same binding performance as the original clone N497/10 but can be detected using standard anti-chicken secondary antibodies allowing flexibility for multiplexing applications. This antibody is expressed recombinantly in HEK293 cells and then affinity purified from the cell culture media and conjugated.</p>
Host	Chicken
Gene ID	55714
Gene Symbol	TENM3
Species	Human, Mouse, Rat
Immunogen	Recombinant protein comprising amino acids aa 161-201 and 343-396 (extracellular and cytoplasmic domains) of mouse ODZ3 (Uniprot#: Q9WTS6) connected by a Gly8 linker, produced recombinantly in E.Coli.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-52434

NBP2-13687PEP	ODZ3 Recombinant Protein Antigen
236-EG-200	EGF [Unconjugated]
6507-IL-010/CF	IL-4 [Unconjugated]
NB110-68136	Tenascin C Antibody (4C8MS) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-52434

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

