

EphB2 Antibody (2D12C6)


Product Description	
Host	Mouse
Gene ID	2048
Gene Symbol	EPHB2
Species	Human
Applications/Dilutions	
Dilutions	ELISA 1:10000 Immunocytochemistry/Immunofluorescence 1:200 - 1:1000 Western Blot 1:500 - 1:2000

P24 well → 100.000 cells/well


POECs → porcine Outgrowth endothelial cells isolated from fetal porcine blood cord

4 days in Endothelial cell medium (Pelotech)


Immunofluorescence → Dil 1:100



pOECs p2



pOECs p2



pOECs p2